



Spartan Chemical Company, Inc., 1110 Spartan Drive, Maumee, OH 43537  
 Toll Free: 1-800-537-8990, Fax: 419-536-8423, Fax-On-Demand: 888-383-4666  
 E-Mail: spartan@spartanchemical.com, Web Site: www.spartanchemical.com  
 ©SCC 2003



re:



**Clothesline**  
**FRESH™**  
**Enzyme Spotter**



### PRODUCT DESCRIPTION:

Clothesline Fresh Enzyme Spotter is a concentrated laundry spotter and wash additive. The unique enzymatic activity of Clothesline Fresh Enzyme Spotter allows for quick penetration, digestion, and removal of the toughest stains and soils. Utilizing Clothesline Fresh Enzyme Spotter as part of a complete laundry care line will reduce rewash and increase linen life. The proprietary blend of enzymes found in Clothesline Fresh Enzyme Spotter allows for the quick removal of many types of stains particularly protein, starch, and manann gum based stains. Examples of protein stains are blood, egg, grass, dairy products, and bodily fluids. Examples of starch based stains include pasta and fruit stains, sugar based stains, sauces, baseball and tennis court red clay. Manann Gums are associative type thickeners and stabilizers that are found in thousands of products. Manann gums provide tackiness that often results in staining when deposited onto clothing, use Clothesline Fresh Enzyme Spotter to attack these staining situations. On most fabrics, Clothesline Fresh Enzyme Spotter will not contribute to fabric degradation.

Use Clothesline Fresh Enzyme Spotter to remove these common spots:

- blood
- egg
- grass
- dairy
- fruit
- red clay
- bodily fluids
- sauces

### COLOR AND NUMBER SYSTEM:

Each Clothesline Fresh laundry product has been color coded and number coded for easy product identification. All product information and literature, plus, each product label, carton and secondary label reflects the individual product's color and number, making training easier and faster.

### STORAGE:

Clothesline Fresh laundry products can be stored in a wide variety of environments and locations. Clothesline Fresh Enzyme Spotter can withstand three freeze/thaw cycles without affecting the performance of the product.

### DIRECTIONS FOR USE:

#### DO NOT MIX WITH OTHER CHEMICALS.

Some individuals may be sensitive to ingredients in this product. Before use, read product label and MSD sheet.

- Always test fabric for color fastness and compatibility in an inconspicuous area.
- Before using a specific spotter, attempt to identify the stain, to ensure efficient stain removal.
- Always wash and rinse fabric thoroughly after spot treatment.

**To Spot Fabric:**

1. Saturate a clean white towel with Clothesline Fresh Enzyme Spotter.
2. Rub stain vigorously with the wet towel.
3. Blot fabric with a separate clean, dry towel.
4. Repeat steps as necessary to achieve desired results.

**As a Pre-Spotter:**

1. Spray Clothesline Fresh Enzyme Spotter liberally onto suspected stains.
2. Allow to remain on stained fabric for at least 5 - 10 minutes before initial laundering.

**As a Pre-Soak:**

Clothesline Fresh Enzyme Spotter can also be used as a presoak for heavily soiled or stained fabric. When used as a Pre-Soak, Clothesline Fresh Enzyme Spotter can be used in concentrated form or diluted down to 1:5 depending on soil load levels.

1. Place stained fabric in a suitable containment vessel.
2. Pour Clothesline Fresh Enzyme Spotter directly on stained fabric or diluted down to 1:5 depending on soil load levels.
3. Allow stained fabrics to soak for at least 15 minutes before initial laundering.

**As an Enzyme Wash Additive:**

Clothesline Fresh Enzyme Spotter can be used as a wash wheel additive to help digest and break down stubborn stains quickly and more efficiently.

1. Add 2 - 3 oz. (¼ to 1/3 cup) per wash load in addition to your regular laundry chemicals.
2. For wash load sizes above 100 lbs., add 4 - 6 oz. of Clothesline Fresh Enzyme Spotter.
3. Extremely heavy stains and soil loads may require an additional wash cycle.

**SPECIFICATION DATA:**

pH – 7.7 – 8.3

Specific Gravity – 1.003 @ 24°C/75°F

Density – 8.347

Stability -

- a. Shelf – One year minimum @ 24°C/75°F
- b. Accelerated – 30 days minimum @ 49°C/120°F
- c. Freeze/Thaw – Can withstand three freeze/thaw cycles

Phosphate free

EDTA free

Nonylphenol ethoxylate free

Silicate free

**PACKAGING:**

Clothesline Fresh Enzyme Spotter is packaged in gallons, four per case and RTU Handi Spray<sup>®</sup> quart bottles, 12 per case with 3 three trigger sprayers. Label copy is available in English, Spanish and French. Secondary labels are also available.

**Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting [www.spartanchemical.com](http://www.spartanchemical.com).**

**GUARANTEE:**

Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of date of manufacture.